



A QIAGEN / BD Company



PAXgene[®]

Bone Marrow RNA Tube

FOR THE COLLECTION OF HUMAN BONE MARROW AND
IMMEDIATE STABILIZATION OF INTRACELLULAR RNA

For Research Use Only

Not for use in diagnostic procedures

To be used only in conjunction with the
PAXgene Bone Marrow RNA Kit (30) (Cat. No. 764133)

I. Intended Use

The PAXgene Bone Marrow RNA System consists of a bone marrow collection tube (PAXgene Bone Marrow RNA Tube) and nucleic acid preparation kit (PAXgene Bone Marrow RNA Kit). It is intended for the collection, storage, and transport of human bone marrow and stabilization of intracellular RNA in a closed tube and subsequent isolation and purification of total RNA.

Product Features

PAXgene Bone Marrow RNA Tube 764114	2.0 ml • 16 x 125 mm Plus • BD Hemogard™	11.1 ml additive 50 tubes/case	Sterile RUO
--	---	-----------------------------------	----------------

For Research Use Only. Not For Use In Diagnostic Procedures. No claim or representation is intended to provide information for the diagnosis, prevention, or treatment of a disease. The performance characteristics of this product have not been fully established.

Do Not Use For Blood Collection.

II. Summary and Explanation

The aspiration and collection of whole bone marrow is the first step in many molecular assays used to study intracellular RNA. A major challenge in this type of testing is the instability of intracellular RNA, which rapidly degrades within hours after bone marrow aspiration. Furthermore, the expression of certain species of RNA, through the process of gene induction, increases in vitro after bone marrow collection. Both in vitro RNA degradation and gene induction can lead to an under- or overestimation of in vivo gene transcript number.

PAXgene Bone Marrow RNA Tubes contain an additive that stabilizes the in vivo gene transcription profile by reducing in vitro RNA degradation and minimizing gene induction. When used in conjunction with the PAXgene Bone Marrow RNA Kit, the PAXgene Bone Marrow RNA Tube provides stabilization and purification of intracellular RNA for accurate detection and quantification of gene transcripts.

III. Specimen Collection and Preparation for Analysis

1. Label PAXgene Bone Marrow RNA Tube appropriately with subject identification.
2. Aspirate bone marrow using standard medical practices. Collect the bone marrow in a syringe with a standard anticoagulant (not supplied).
3. Remove the needle from the syringe used to collect bone marrow. Remove the BD Hemogard Closure from the PAXgene Bone Marrow RNA Tube as stated in Section VII. Add 2.0 ml bone marrow to the labeled PAXgene Bone Marrow RNA Tube as soon as possible after bone marrow aspiration. Use the gradations on the syringe to measure the amount of bone marrow. Alternatively, transfer the bone marrow to a secondary tube and pipet 2.0 ml bone marrow into the PAXgene Bone Marrow RNA Tube.

Note: Use of greater than 2.0 ml will interfere with the stabilization properties of the PAXgene Bone Marrow RNA Tube, leading to reduced RNA yields and degraded RNA. Use of less than 2.0 ml can lead to incorrect analytical results or poor product performance.

Note: The bone marrow sample is subject to RNA degradation and changes in gene expression until the bone marrow is transferred into the PAXgene Bone Marrow RNA Tube. To reduce changes of transcript levels that occur within minutes of aspiration, reduce as much as possible the time between bone marrow aspiration and transfer to the PAXgene Bone Marrow RNA Tube.

4. Close the PAXgene Bone Marrow RNA Tube with the original BD Hemogard Closure as stated in Section VIII. To close the tube, replace closure over tube. Next, twist and push down firmly until stopper is fully reseated. Complete reinsertion of the stopper is necessary for the closure to remain securely on the tube and prevent leakage during handling.
5. Immediately after bone marrow transfer and recapping, secure closure with thumb while holding tube in one hand, and mix the tube contents vigorously by powerfully inverting the PAXgene Bone Marrow RNA Tube 8–10 times.

Note: For sufficient RNA stabilization, it is important to mix bone marrow and tube solution completely. Some anticoagulants (e.g. lithium heparin) are known to complicate efficient mixing if tubes are only gently inverted.

6. Store the PAXgene Bone Marrow RNA Tube upright at room temperature (18–25°C) for a minimum of 2 hours before processing or transferring to refrigerator (2–8°C) or freezer (–20°C). Bone Marrow samples filled in PAXgene Bone Marrow RNA Tubes are stable at room temperature (18–25°C), refrigerated (2–8°C), or frozen (–20°C / –70°C).
 - Stand the PAXgene Bone Marrow RNA Tubes upright in a wire rack. Do not use a polystyrene foam tray as this may cause tubes to crack during freezing.
 - If PAXgene Bone Marrow RNA Tubes are to be stored at temperatures lower than –20°C, freeze them first for 24 hours at –20°C, then transfer them to –70°C.
 - Frozen PAXgene Bone Marrow RNA Tubes are subject to breakage upon impact. To reduce the risk of breakage during transport, frozen tubes should be treated in the same manner as glass tubes. Users must validate their own freezing and shipping protocol for PAXgene Bone Marrow RNA Tubes.

Note: Currently available data shows stabilization of intracellular RNA for 13 months at –20°C. For more information from ongoing studies evaluating stability for longer time periods or at other temperatures, please contact QIAGEN Technical Services (techservice-eu@qiagen.com).

Procedure for Specimen Preparation for Analysis

Bone marrow samples must be processed according to the instructions provided with the PAXgene Bone Marrow RNA Kit (Cat. No. 764133).

IV. Limitations of System

1. The PAXgene Bone Marrow RNA System (Tube and RNA Preparation Kit) is for research use only, the use of this product for diagnostic procedures and patient management is strictly prohibited. Neither the clinical utility nor the performance of the PAXgene Bone Marrow RNA Tube and Kit as part of an in vitro diagnostic procedure have been established.
2. Bone marrow samples obtained with the PAXgene Bone Marrow RNA Tube should be prepared only with the PAXgene Bone Marrow RNA Kit (Cat. No. 764133).
3. Over-filling or under-filling of tubes will result in an incorrect bone marrow-to-additive ratio and may lead to incorrect analytic results or poor product performance.
4. Endotoxin not controlled. Bone marrow and bone marrow components collected and processed in the tube are not intended for infusion or introduction into the human body.


V. Warnings and Precautions

1. Contents of this tube may cause irritation to eyes, respiratory system and skin.
 - After inhalation, supply fresh air, and seek medical advice in case of complaints.
 - After skin contact, immediately wash with water and soap and rinse thoroughly.
 - After eye contact, rinse immediately with plenty of water for at least 15 minutes and seek medical advice.
 - After swallowing, immediately seek medical advice.
2. Practice Universal Precautions. Use gloves, gowns, eye protection, other personal protective equipment, and engineering controls to protect from bone marrow splatter, bone marrow leakage, and potential exposure to blood borne pathogens.
3. Handle all biologic samples and bone marrow aspiration devices (e.g., syringes or other collection devices) according to the policies and procedures of your facility. Obtain appropriate medical attention in the event of any exposure to biologic samples (for example, through a puncture injury), since they may transmit viral hepatitis, HIV, or other infectious diseases. Utilize any built-in used-needle protector, if the aspiration device provides one. PreAnalytiX does not recommend re-shielding used needles; however, the policies and procedures of your facility may differ and must always be followed.
4. Discard all bone marrow collection tubes and accessories in biohazard containers approved for disposal.
5. Do not draw blood using the PAXgene Bone Marrow RNA Tube.
6. Do not reuse PAXgene Bone Marrow RNA Tubes.
7. Do not transfer a sample from a syringe to a tube through the stopper.
8. Excessive centrifugation speed (over 10,000 RCF) may cause PAXgene Bone Marrow RNA Tube breakage, exposure to bone marrow and possible injury.

VI. Storage

Store the unused PAXgene Bone Marrow RNA Tubes at 4–25°C.

VII. Instructions for Removal of BD Hemogard Closure




1. Grasp the PAXgene Bone Marrow RNA Tube with one hand, placing the thumb under the BD Hemogard Closure (for added stability, place arm on solid

surface). With the other hand, twist the BD Hemogard Closure while simultaneously pushing up with the thumb of the other hand **ONLY UNTIL THE TUBE STOPPER IS LOOSENED**.

2. Move thumb away before lifting closure. **DO NOT** use thumb to push closure off tube. Caution: If the tube contains blood, an exposure hazard exists. To help prevent injury during closure removal, it is important that the thumb used to push upward on the closure is removed from contact with the tube as soon as the BD Hemogard Closure is loosened.
3. Lift closure off tube. In the unlikely event of the plastic shield separating from the rubber stopper, **DO NOT REASSEMBLE CLOSURE**. Carefully remove rubber stopper from tube.

VIII. Instructions for Reinsertion of BD Hemogard Closure



1. Replace closure over tube.
2. Twist and push down firmly until stopper is fully resealed. Complete reinsertion of the stopper is necessary for the closure to remain securely on the tube during handling.

IX. Technical Assistance

If you have any questions regarding the PAXgene Bone Marrow RNA System, contact one of the QIAGEN Technical Services Departments listed in the QIAGEN – Customer Service section.

X. Ordering Information

PAXgene Products

PAXgene Bone Marrow RNA Tubes (50) / 764114

Contents: 50 Bone Marrow Collection Tubes.

Use in conjunction with the PAXgene Bone Marrow RNA Kit (30).

PAXgene Bone Marrow RNA Kit (30) / 764133

Content: 30 PAXgene RNA Spin Columns, 30 PAXgene Shredder Spin Columns, Processing Tubes, RNase-free DNase-I, RNase-free Reagents and Buffers.

Use in conjunction with the PAXgene Bone Marrow RNA Tubes (50).

QIAGEN – Customer Service

Australia

Orders: 1.800.243.800
Technical: 1.800.243.066
Fax: 03.9840.9888

Austria

Orders: 0800.28.10.10
Technical: 0800.28.10.11
Fax: 0800.28.10.19

Belgium

Orders: 0800.79612
Technical: 0800.79556
Fax: 0800.79611

Canada

Orders: 800.572.9613
Technical: 800.DNA.PREP (800.362.7737)
Fax: 800.713.5951

China

Orders: 86.21.3865.3865
Technical: 800.988.0325
Fax: 86.21.3865.3965

Denmark

Orders: 80.885945
Technical: 80.885942
Fax: 80.885944

Finland

Orders: 0800.914416
Technical: 0800.914413
Fax: 0800.914415

France

Orders: 01.60.920.926
Technical: 01.60.920.930
Fax: 01.60.920.925

Germany

Orders: 02103.29.12000
Technical: 02103.29.12400
Fax: 02103.29.22000

Hong Kong

Orders: 800.933.965
Technical: 800.930.425
Fax: 800.930.439

Ireland

Orders: 1800.555.049
Technical: 1800.555.061
Fax: 1800.555.048

Italy

Orders: 800.789.544
Technical: 800.787.980
Fax: 02.334304.826

Japan

Orders: 03.6890.7300
Technical: 03.6890.7300
Fax: 03.5547.0818

Korea (South)

Orders: 080.000.7146
Technical: 080.000.7145
Fax: 02.2626.5703

Luxembourg

Orders: 8002.2076
Technical: 8002.2067
Fax: 8002.2073

Mexico

Orders: 01.800.7742.639
Technical: 01.800.7742.436
Fax: 01.800.1122.330

The Netherlands

Orders: 0800.0229592
Technical: 0800.0229602
Fax: 0800.0229593

Norway

Orders: 800.18859
Technical: 800.18712
Fax: 800.18817

Singapore

Orders: 1800.742.4362
Technical: 1800.742.4368
Fax: 65.6854.8184

Spain

Orders: 91.630.7050
Technical: 91.630.7050
Fax: 91.630.5145

Sweden

Orders: 020.790282
Technical: 020.798328
Fax: 020.790582

Switzerland

Orders: 055.254.22.11
Technical: 055.254.22.12
Fax: 055.254.22.13

BD - Customer Service**Argentina, Uruguay and Paraguay**

Orders: 0800.444.5523
E-mail: crc_argentina@bd.com

Australia

Orders: 1.800.656.100
Fax: 1.800.656.110
E-mail: bd_anz@bd.com

Austria

Orders: 43.1.706.36.60
Fax: 43.1.706.36.60-11
Email: customerservice.vac.ga@europe.bd.com

Belgium

Orders: 32.53.720.556
Fax: 32.53.720.549
Orders: orders.be@europe.bd.com

Brazil

Orders: 0800.055.56.54
E-mail: consultoria_vacutainer@bd.com

Canada

Orders: 800.268.5430
Fax: 800.565.0897
E-mail: custsvccan@bd.com

Central and Eastern Europe (CEE)

Orders: 48.22.377.11.11
Fax: 48.22.377.11.02
Bulgaria: info_bulgaria@europe.bd.com
Czech Republic: info_czech@europe.bd.com
Croatia: info_croatia@europe.bd.com
Hungary: info_hungary@europe.bd.com
Poland: info_poland@europe.bd.com
Romania: info_romania@europe.bd.com
SEE: info_balkan@europe.bd.com
Serbia: info_serbia@europe.bd.com

UK

Orders: 01293.422.911
Technical: 01293.422.999
Fax: 01293.422.922

USA

Orders: 800.426.8157
Technical: 800.DNA.PREP (800.362.7737)
Fax: 800.718.2056

Slovakia: info_slovakia@europe.bd.com
Slovenia: info_slovenia@europe.bd.com
Technical: bd_service_poland@bd.com

Denmark

Orders: 45.43.43.45.66
Fax: 45.43.96.56.76
Orders: ordre.dk@europe.bd.com
Technical: bddenmark@europe.bd.com

Finland

Orders: 358.9.88.70.780
Fax: 358.9.88.70.7816
Orders: tilaukset.fi@europe.bd.com
Technical: bdsuomi@europe.bd.com

France

Orders: 33.476.68.94.96
Fax: 33.476.68.36.93
Orders: commandesfr@europe.bd.com
E-mail: serviceclientbdf@europe.bd.com
Technical: vacutainerfr@bd.com

Germany

Orders: 49.6221.305.553
Fax: 49.6221.305.379
Email: customerservice.vac.ga@europe.bd.com

Italy

Orders: 39.02.48.240.500
Technical: 39.335.424.388
Fax: 39.02.40.91.80.11
E-mail: ordini.it@europe.bd.com

Middle East & Africa (EMA)

Orders: 971.4.3379525
Fax: 971.4.03379551
E-mail: EMA_PAS@bd.com

The Netherlands

Orders: 31.20.582.94.20

Fax: 31.20.582.94.21

Orders: orders.nl@europe.bd.com

New Zealand

Orders: 0800.572.468

Fax: 0800.572.469

E-mail: nz_customerservice@bd.com

Norway (Puls Handicare)

Orders: 47.23.32.30.00

Fax: 47.23.32.30.99

E-mail: kundeservice@puls-norge.no

Spain, Portugal and Andorra

Orders: 34.91.848.8104

Technical: 34.902.27.17.27

Fax: 34.91.828.81.48

E-mail: info.spain@europe.bd.com

Sweden

Orders: 46.8.775.51.60

Fax: 46.8.775.51.95

Orders: order.se@europe.bd.com

Technical: bds sweden@europe.bd.com

Switzerland

Orders: 41.61.485.22.22

Fax: 41.61.485.22.00

E-mail: infoch@europe.bd.com

UK

Orders: 44.1.865.781.666

Fax: 44.1.865.717.313

E-mail: bduk_customerservice@europe.bd.com

USA

Orders: 888.237.2762

Technical: 800.631.0174

Fax: 800.847.2220

Technical: www.bd.com/vacutainer/contact

Symbol and Mark Key



Batch Code



Catalog Number



Contains sufficient for <n> tests



Keep Away from Sunlight



Manufacturer



Method of Sterilization Using
Irradiation



Temperature Limitation



Recyclable



Use By



H315 – Causes skin irritation.

H317 – May cause an allergic skin reaction.




H318 – Causes serious eye damage.



H411 – Toxic to aquatic life with long last effects.

The PAXgene Bone Marrow RNA Tube is made in the UK by BD for PreAnalytiX GmbH.

PAXgene and PreAnalytiX are trademarks of PreAnalytiX GmbH; QIAGEN and QIAamp are trademarks of QIAGEN NV; all other trademarks are property of Becton, Dickinson and Company.

 PreAnalytiX GmbH, Feldbachstrasse, 8634 Hombrechtikon, CH ©2016 PreAnalytiX GmbH