Influence of Centrifugation Conditions for Plasma Processing on ccfDNA Yield ¹Daniel Groelz, ¹Tomasz Krenz, ¹Eric Provencher and ¹Thorsten Voss

¹PreAnalytiX GmbH, Hombrechtikon, Switzerland


Introduction

Circulating, cell-free DNA (ccfDNA) in blood has become widely accepted for biomarker analysis in academic and clinical research and beyond. Considering the usual low sample input and required high sensitivity for target detection, optimal outcomes require validation and control of the entire workflow. The PAXgene Blood ccfDNA System (RUO) integrates preanalytical steps from sample collection to ccfDNA extraction for subsequent analysis by quantitative PCR (qPCR), digital PCR (ddPCR) and next-generation sequencing (NGS).

Here we present the influence of centrifugation conditions for plasma processing on ccfDNA yield and compatibility with primary tube handling on the QIAsymphony[®] SP instrument using the QIAsymphony PAXgene[®] Blood ccfDNA Kit* and corresponding protocols.

*For Research Use Only. Not for use in diagnostic procedures.

Stability claims for the PAXgene Blood ccfDNA Tube are 37°C for up to 1 day and 2–30°C for up to 7 days.



. Option A is the workflow described in the QIAsymphony PAXgene Blood ccfDNA Kit Handbook using the optional second centrifugation step. Option B is the new primary tube handling workflow using the PAXgene Blood ccfDNA Tube to be processed directly on the QIAsymphony SP instrument. Time savings based on processing a 24 sample batch.

Methods


Blood from healthy consented donors was drawn into PAXgene Blood ccfDNA Tubes (PreAnalytiX®)* or Cell-Free DNA BCT[®] (Streck[®]) and was processed into plasma directly after phlebotomy or stored for 6–7 days at room temperature (15–25°C) to simulate a typical processing delay in a routine setting. Different plasma processing protocols were tested in 4 studies (Table 1).

The level of hemolysis in the separated plasma was analyzed by photometric measurement of the absorption at 414 nm, which is the absorption maximum for hemoglobin. ccfDNA was isolated from PAXgene plasma on the QIAsymphony SP (QIAGEN[®]) automated platform using the corresponding QIAsymphony PAXgene Blood ccfDNA Kit* and protocol and from Streck plasma using the QIAsymphony DSP Circulating DNA Kit (QIAGEN) and protocol. Yield and stability of ccfDNA generated in Studies 1–3 were measured by qPCR (QIAGEN Rotor-Gene® Q) and ddPCR (BioRad QX 200 DR Droplet Digital PCR System) with assays amplifying a 66 bp fragment of the 18S rDNA sequence. Compatibility with NGS was evaluated by sequencing ccfDNA samples from Study 4 on the QIAGEN GeneReader™ NGS System, including the GeneRead™ QIAact Actionable Insights Tumor (AIT) Panel* for PCR target enrichment and library preparation on the QIAcube[®] instrument, QC with capillary electrophoresis, sequencing on the GeneReader instrument* and data management with the QIAGEN Clinical Insight (QCI[™]) Analyze* tool.

*For Research Use Only. Not for use in diagnostic procedures.


Study 1 Brake settings		25 donors	
	Blood collection & storage	10 PAXgene Blood ccfDNA Tubes	
		5 tubes: direct processing (t0), 5 tubes: 6 days at 25°C	
	Plasma processing (Sigma 4K15C centrifuge)	Two centrifugation steps	15´at 1,900 × g, 10´ at 1,900 × g Maximum brake (19) 15´ at 1,900 × g, 10´ at 1,900 × g Medium brake (9)
		One centrifugation step	15´ at 1,900 × g Maximum brake (19)
			15´ at 1,900 × g Medium brake (9)
			15´ at 1,900 × g No brake (0)
	ccfDNA extraction	4.8 ml plasma	QIAsymphony PAXcircDNA_STA_4800
Study 2 Centrifugation forces & times	Blood collection & storage	25 donors	
		8 PAXgene Blood ccfDNA Tubes, 2 Streck Cell-Free DNA glass tubes	
		4 PAXgene tubes: direct processing (t0), 4 PAXgene & 2 Streck tubes: 7 days at 25°C
	PAXgene plasma processing (Sigma 4K15C centrifuge)	Two centrifugation steps	15´ at 1,900 × g, 10´ at 1,900 × g Medium brake (9)
		One centrifugation step	15´ at 1,900 × g Medium brake (9)
			10´ at 2,000 × g Medium brake (9)
			5´ at 3,000 × g Medium brake (9)
	Streck cfDNA spin 1 & spin 2 protocols* (Sigma 4K15C centrifuge)	Two centrifugation steps	Spin 1: 20´ at 300 × g, 10´ at 5,000 × g Medium brake (9)
			Spin 2: 10´ at 1,600 × g, 10´ at 16,000 × g Medium brake (9)
	ccfDNA extraction	PAXgene plasma: 2.4 ml	QIAsymphony PAXcircDNA_STA_2400
		Streck tubes: 2.0 ml	QIAsymphony circDNA_2000
Study 3 Primary tube work- flow	Blood collection & storage	17 donors	
		10 PAXgene Blood ccfDNA Tubes	
		5 tubes: direct processing (t0), 5 tubes: 7 days at 25°C	
	Plasma processing (Sigma 4K15C centrifuge)	One centrifugation Medium brake (9)	15´ at 1,900 × g
			5´ at 1,600 × g
			15´ at 1,600 × g
			5´ at 3,000 × g
			15´ at 3,000 × g
	ccfDNA extraction	4.0 ml / 2.4 ml	QIAsymphony primary tube PAXcircDNA_STA_4000 or 2400
Study 4 NGS workflow	Blood collection & storage	60 donors	
		10 PAXgene Blood ccfDNA Tubes	
		1× direct processing, 1× 7 days at 25°C	
	Plasma processing (Sigma 4K15C centrifuge)	Two centrifugation steps	15´ at 1,900 × g, 10´ at 1,900 × g Medium brake (9)
	ccfDNA extraction	4.8 ml plasma	QIAsymphony PAXcircDNA_STA_4800

*Protocols according to IFU for Streck Cell-Free DNA BCT




For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN and PreAnalytiX kit handbook or user manual. Handbooks and user manuals are available at www.qiagen.com and www.preanalytix.com or can be requested from QIAGEN Technical Services or your local distributor. Trademarks are the property of their respective owners. © 2018 PreAnalytiX GmbH. PreAnalytiX, the PreAnalytiX Logo and all other trademarks are property of PreAnalytiX GmbH, Hombrechtikon, CH.

Results



- blood draw.
- 25 (>80%).
- identified in healthy donors using the AIT panel.



- Centrifugation parameters are important for optimal plasma quantity and quality and are dependent on the tube type used for blood collection
- Main findings from the different centrifugation conditions tested:
- High brake settings can lead to angled plasma/cell border Low centrifugation forces and times (under 1,600 × g for less than 15 min) can lead to reduced plasma yield and insufficient separation of plasma from cellular fraction, which can lead to gDNA contamination of ccfDNA by transferred white blood cells
- No tube breakage was observed for PAXgene Blood ccfDNA Tubes Centrifugation in range of 1,600 × g to 3,000 × g for 15 min has no impact on hemolysis and results in no significant differences for total yield and in situ stability (yield change over storage time)



Study 2 – Influence of different centrifugation conditions during plasma processing

- Recommended protocol:
- 15 minutes at 1,600–3,000 × g with medium brake setting - Optional: Second centrifugation, 10 minutes at 1,600–3,000 × g with medium brake setting
- Data from the first 3 studies imply that the recommended centrifugation parameters compared to the procedure described in the product literature will not impact:
- Plasma yield and quality
- Compatibility with automated ccfDNA isolation procedures including primary tube handling workflow
- Downstream performance with qPCR and digital droplet PCR Eluates from the PAXgene Blood ccfDNA System are compatible with the GeneReader NGS workflow