Verification of a Complete Sample to Insight[®] Liquid Biopsy Workflow – NGS of ccfDNA From Stabilized Blood Tomasz Krenz, Andrea Ullius, Ricardo Huebel, Thorsten Voss, Eric Provencher, Timothy R. Buirkle and Daniel Groelz

PreAnalytiX GmbH, Hombrechtikon, Switzerland

Introduction

Analysis of circulating cell-free DNA (ccfDNA) from blood in the fields of life science, clinical research and beyond is expanding and has become widely accepted. The complete workflow entails collection and stabilization of blood, transport, plasma processing, extraction of ccfDNA and downstream analysis. Considering the usual low sample input and required high sensitivity for quality control (QC), optimal outcomes require verification of the entire workflow.

In this verification study, human blood was collected, stabilized, and then processed to isolate ccfDNA from plasma using the PAXgene® Blood ccfDNA System for subsequent next-generation sequencing (NGS) on the QIAGEN GeneReader[®] system. To verify the workflow, several sample QC criteria were assessed at each step of the workflow, including ccfDNA yield, target enrichment, library preparation, target sequencing and analysis.

Methods

Blood from 60 healthy consented donors was collected into PAXgene Blood ccfDNA Tubes* and stored for 7 days at 25°C to simulate a typical processing delay in a routine setting. Automated ccfDNA extraction was performed on the QIAsymphony[®] SP instrument using the QIAsymphony PAXgene Blood ccfDNA Kit* and protocol. ccfDNA stability after blood storage was confirmed by qPCR. ccfDNA samples were sequenced on the QIAGEN GeneReader NGS System, including the GeneRead™ QIAact Actionable Insights Tumor (AIT) Panel* for PCR target enrichment, library preparation on the QIAcube[®] instrument, QC with capillary electrophoresis, sequencing on the GeneReader instrument* and data management with the QIAGEN Clinical Insight (QCI™) Analyze* tool.

Target Enrichment & Library Preparation



• All samples (60/60) passed the required QC criteria after target enrichment (amplicon size <170 bp) and library preparation (increase of amplicon size by 92 bp (adapters) and absence of unspecific products).

Compatibility of the liquid biopsy workflow in target enrichment (A) and library preparation (B) NGS procedures. Blood from 60 healthy donors was stored in PAXgene Blood ccfDNA Tubes for 7 days at 25°C prior to ccfDNA extraction. Maximum volume of ccfDNA eluate was used for the upfront NGS workflow using the GeneReader workflow and the GeneRead QIAact Actionable Insights Tumor (AIT) panel. The products after target enrichment (A) and library preparation (B) were analyzed for correct fragment size (I.) and total product concentration (II.) using capillary electrophoresis on the QIAxcel instrument. Representative electropherograms and single data points are shown. The red line indicates the minimal DNA concentration requirement for the next NGS step.

- Blood from 60 healthy donors



- quality 25 (>80%).



Study Design – ccfDNA NGS Workflow

• PAXgene Blood ccfDNA Tube, QIAsymphony PAXgene Blood ccfDNA Kit, GeneReader NGS system: Sample to Insight



AACR Annual Meeting 2018 Poster number 459

In Situ ccfDNA Stability

• After blood storage for 7 days, ccfDNA yield was comparable to ccfDNA yield observed directly after blood draw.