



How to Collect Blood Using the PAXgene® Blood ccfDNA Tube (RUO)

For Circulating Cell-Free DNA (ccfDNA)




Before Blood Collection:

- 1a. Ensure that the PAXgene Blood ccfDNA Tube* (A) is at room temperature (15–25°C) prior to use and is properly labeled with specimen identification.
- 1b. The PAXgene Blood ccfDNA Tube should be the last tube drawn in the phlebotomy procedure.




Venipuncture:

2. Using a blood collection set such as the BD Vacutainer UltraTouch Push Button or BD Vacutainer Safety-Lok Blood Collection Set (B) and the BD Vacutainer One-Use Holder (C), perform venipuncture using your institution's recommended procedure for standard venipuncture.



Blood Collection:

- 3a. Hold the PAXgene Blood ccfDNA Tube vertically below the blood donor's arm during blood collection.
- 3b. Ensure that the blood has stopped flowing into the tube before removing the tube from the holder.



After Blood Collection:

- 4a. Gently invert the PAXgene Blood ccfDNA Tube 8 times.
- 4b. Store the PAXgene Blood ccfDNA Tube upright until centrifugation and plasma processing.† Blood-filled PAXgene Blood ccfDNA Tubes can be stored for up to 10 days at temperatures up to 25°C, 7 days at temperatures up to 30°C, or 3 days at temperatures up to 37°C. Do not store blood-filled tubes below 2°C.

Before Blood Collection



After Blood Collection



Ref #	Product Description
768115	PAXgene Blood ccfDNA Tube (RUO)
367364	BD Vacutainer® UltraTouch™ Push Button Blood Collection Set
367281 (US) 367286 (CE)	BD Vacutainer Safety-Lok™ Blood Collection Set
364815	BD Vacutainer One-Use Holder

Visit www.preanalytix.com for more information.

* For Research Use Only
Not for use in diagnostic procedures

www.PreAnalytiX.com

† Refer to PAXgene Blood ccfDNA Tube Handbook.
PAXgene and PreAnalytiX are registered trademarks of PreAnalytiX, Hombrechtikon, CH
All other trademarks are property of Becton, Dickinson and Company.
© 2020 PreAnalytiX GmbH

PreAnalytiX

A QIAGEN / BD Company