

Preanalytical workflow enabling cfDNA analysis from urine samples

¹Daniela Mancarella and ¹Thorsten Voss
¹PreAnalytiX GmbH, Hombrechtikon, Switzerland
✉ Please contact Daniela.Mancarella@qiagen.com for further questions or go to www.preanalytix.com for information on PreAnalytiX


Introduction

Urine has become an important source of information in the liquid biopsy field. In contrast to blood, specifications for the collection, storage, transport and processing of urine intended for molecular examination are not widely established. Preanalytical specifications were published for urine cell-free DNA (cfDNA) only recently (CEN/TS 17811:2022). In this study, we investigated post-collection changes to cfDNA profiles in urine samples and present the performance of an optimized preanalytical workflow for cfDNA analysis.

Methods

Urine from apparently healthy, consented female individuals was collected. Where indicated in the figures, urine specimens were pooled (urine of at least 3 individuals). Urine was spiked with cell-free supernatant from male urine or with a cfDNA standard containing *PIK3CA* mutations (Multiplex 5% allele frequency (AF) cfDNA standard, SensiD). The urine was stabilized with a urine collection/stabilization solution under development at PreAnalytiX or was left unstabilized. Urine samples were stored for varying durations at different temperatures in time course experiments to simulate research settings. After storage, urine samples were centrifuged to remove cells and cfDNA was isolated from the supernatant via the QIAAsymphony® platform (QIAGEN). Autosomal (cfDNA target 1) and male-specific targets (cfDNA target 2) were quantified using the Rotor-Gene® Q instrument (QIAGEN) and the Investigator® Quantiplex® Pro RGQ assay (QIAGEN). *PIK3CA* H1047R mutation was analyzed by a *PIK3CA* qPCR assay on the Rotor-Gene Q instrument or by dPCR on the QIAcuity® Digital PCR System using the dPCR LNA *PIK3CA* Mutation Assay (QIAGEN). Fragment size distribution was determined by TapeStation® Cell-free DNA ScreenTape® (Agilent Technologies).



Experimental setup for cfDNA analysis using stabilized urine





Results

Storage-induced cfDNA profile changes can be prevented by urine stabilization

Degradation of cfDNA and release of genomic DNA is minimized by urine stabilization




(B) Time course of cfDNA spike-in in unstabilized and stabilized urine over 10 days of storage



Preservation of qPCR- and dPCR- detectable cfDNA targets

cfDNA from urine stored after stabilization can be detected with qPCR and dPCR



Conclusion

This study showed that post-collection changes in unstabilized urine result in an artificial cfDNA profile and loss of target of interest leading to failed outcomes. Urine stabilization with a collection/stabilization solution under development at PreAnalytiX minimized DNA degradation and gDNA release. It enables urine storage, allowed analysis of urine cfDNA profile and target detection by qPCR and dPCR.

Disclaimer

Product names mentioned are MBA and are for Molecular Biology Applications. They are not for use in diagnosis, prevention, or treatment of a disease. Some products may not be available in all countries. For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN or PreAnalytiX handbook or user manual. QIAGEN and PreAnalytiX handbooks and user manuals are available at www.qiagen.com or www.preanalytix.com or can be requested from QIAGEN Technical Services or your local distributor. Trademarks are the property of their respective owners. Trademarks: QIAGEN®, QIAcuity®, QIAAsymphony®, Investigator®, Quantiplex®, Rotor-Gene® (QIAGEN Group); ScreenTape®, TapeStation® (Agilent Technologies) © 2024 PreAnalytiX GmbH. PreAnalytiX, the PreAnalytiX Logo and all other trademarks are property of PreAnalytiX GmbH, Hombrechtikon, CH.