

QIASymphony PAXgene[®] Blood RNA Kit Now Available

New kit allows automated purification of intracellular RNA from stabilized blood

PreAnalytiX GmbH, a QIAGEN/BD company, has launched the QIASymphony PAXgene Blood RNA Kit.

The QIASymphony PAXgene Blood RNA Kit provides fully automated purification of intracellular RNA on the QIASymphony SP. Used with PAXgene Blood RNA Tubes (available from BD and BD authorized distributors, cat. no. 762165) for sample collection and stabilization, the system provides a complete preanalytical solution for sample collection, stabilization and purification of high quality RNA using magnetic particle technology. Pre-filled cartridges minimize hands-on time, and bar code labeling enables complete reagent tracking.

For more information on the QIASymphony PAXgene Blood RNA Kit click here:

<http://www1.qiagen.com/Products/QIASymphonyPAXgeneBloodRNAKit.aspx>

The QIASymphony PAXgene Blood RNA Kit is for research use only. No claim or representation is intended to provide information for the diagnosis, prevention or treatment of a disease.